

**Reagents for the EP assays for Unfractionated Heparin**

Anti IIa	
Range 0,005 -0,03 IU/ml	
<b>Buffer R1</b>	
<b>5D Tris-NaCL-EDTA-PEG Buffer pH 8.4</b>	Ref. HB0187
Package	Pouch
Reconstitution	dissolve pouch in water and make up to 1000 ml
Concentration	0.050 M Tris buffer pH 8.4 at 25°C 0.175 M NaCl 0.0075 M EDTA 0.10% (w/v) PEG-6000
<b>Human Thrombin Reagent R2</b>	
<b>5D Human Thrombin</b>	Ref. HE0160
Package	100 IU per vial
Reconstitution	2 ml H2O
Stock concentration	50 IU/ml
Dilute appropriate volume of stock 1:10 with Buffer R1	
<b>Final Concentration</b>	<b>5 IU/ml</b>
<b>Antithrombin III solution R5</b>	
<b>5-D Human Antithrombin</b>	Ref. HE0162
Package	10 IU per vial
Reconstitution	2 ml H2O
Stock concentration	5 IU/ml
Dilute appropriate volume of stock 1:40 with Buffer R1	
<b>Final Concentration</b>	<b>0,125 IU/ml</b>
<b>Chromogenic Substrate for factor IIa</b>	
<b>5D - Chromogenic Substrate for Thrombin</b>	Ref. HS0170
structure:	D-Phe-Pip-Arg-pNA
Package	25 mg per vial 40 µmol/vial (approximatly)
Reconstitution	8 ml H2O
Stock concentration	5 mM
Dilute appropriate volume of stock 1:4 with H2O	
<b>Concentration</b>	<b>1,25 mM</b>

Anti Xa	
Range 0.03 -0,375 IU/ml	
<b>Buffer R1</b>	
<b>5D Tris-NaCL-EDTA-PEG Buffer pH 8.4</b>	Ref. HB0187
Package	Pouch
Reconstitution	dissolve pouch in water and make up to 1000 ml
Concentration	0.050 M Tris buffer pH 8.4 at 25°C 0.175 M NaCl 0.0075 M EDTA 0.10% (w/v) PEG-6000
<b>Bovine FXa Reagent R2</b>	
<b>5D Bovine Xa</b>	Ref. HE0161
Package	30 µg per vial
Reconstitution	2 ml H2O
Stock concentration	15 µg/ml
Dilute appropriate volume of stock 1:5 with Buffer R1	
<b>Final Concentration</b>	<b>3 µg/ml</b>
(Absorbance blanc in 1 cm cuvette at 405 nm between 0.65 and 1.25)	
<b>Antithrombin III solution R6</b>	
<b>5-D Human Antithrombin</b>	Ref. HE0162
Package	10 IU per vial
Reconstitution	2 ml H2O
Stock concentration	5 IU/ml
Dilute appropriate volume of stock 1:5 with Buffer R1	
<b>Final Concentration</b>	<b>1 IU/ml</b>
<b>Chromogenic Substrate for factor Xa</b>	
<b>5D - Chromogenic Substrate for Factor Xa</b>	Ref. HS0171
structure:	D-Arg-Gly-Arg-pNA
Package	25 mg per vial 35 µmol/vial (approximatly)
Reconstitution	7 ml H2O
Stock concentration	5 mM
Dilute appropriate volume of stock 1:5 with H2O	
<b>Concentration</b>	<b>1 mM</b>