Instructions for Use

MRX Specialty Normal Control, K5021

Contents

| | |
|---------------|----|
| English [EN] | 3 |
| Dansk [DA] | 6 |
| Deutsch [DE] | 9 |
| Español [ES] | 12 |
| Suomi [FI] | 15 |
| Magyar [HU] | 18 |
| Italiano [IT] | 21 |
| Norsk [NO] | 24 |
| Svenska [SV] | 27 |

Instructions for Use [EN]

MRX Specialty Normal Control

REF K5021

For *In vitro* Diagnostic Use.

1 Intended use

For quality control of the assays MRX Red Free Protein S, MRX Blue Free Protein S, MRX Protein C and MRX Green vWF Antigen. Intended to be used by professional laboratory personnel using coagulation analysers.

2 Background and principle of method

For background information and principle of the method, refer to the Instructions for Use associated with each individual assay (K5009, K5010, K5020, K5052).

3 Components

MRX Specialty Normal Control consists of:
10 × 1 mL lyophilised citrated human plasma.

4 Metrological traceability

Each new lot of MRX Specialty Normal Control is assigned with each corresponding MRX assay against calibrators with traceability to WHO International Standards, NIBSC code: 03/228 (Free Protein S), 02/342 (Protein C) and 07/316 (vWF Antigen). Refer to the Certificate of Analysis for the lot-specific concentration of each individual analyte.

Specifications for each analyte are listed below:

| Analyte | Specification |
|-----------------|---------------|
| Free Protein S* | 80–110% |
| Protein C | 85–110% |
| vWF Antigen | 90–140% |

*MRX Blue Free Protein S is used for Free Protein S assignment.

5 Warnings and precautions

Wear suitable clothing for protection. Avoid contact with skin and eyes. Do not empty into drains. Waste must be disposed of in accordance with local regulations.

The control contains material of human origin. Each donor has been tested by approved methods and found negative for the presence of HBsAg and anti-HIV I & II and anti-HCV. However, as no method can offer complete assurance

that infectious agents are absent, this material should be handled as potentially infectious.

6 Preparation

- Before opening, carefully tap the vial against a surface to collect the lyophilised material at the bottom.
- Add 1.00 mL deionised water (e.g. MRX Laboratory Water, K5036). The water temperature should be 15 - 25 °C.
- Reseal the vial and let it stand for approximately 15 minutes at 15 - 25 °C.
- Gently mix by swirling or rotating until the content is completely reconstituted.

7 Storage and stability

Store at 2 - 8 °C. After reconstitution, stable for 24 hours at 2 - 25 °C in the closed original vial, provided no contamination occurs.

8 Specimen collection and preparation

For specimen collection and preparation, refer to the Instructions for Use associated with each individual assay (K5009, K5010, K5020, K5052).

9 Procedure

For description of procedure, refer to the Instructions for Use associated with each individual assay (K5009, K5010, K5020, K5052).

10 Material required but not provided

Coagulation analyser, pipettes and the following:

| Reagent | REF |
|-------------------------|-------|
| MRX Red Free Protein S | K5009 |
| MRX Blue Free Protein S | K5010 |
| MRX Protein C | K5020 |
| MRX Green vWF Antigen | K5052 |

| | |
|---|-------|
| Control material | REF |
| MRX Specialty Abnormal Control | K5022 |
| Solutions | REF |
| Deionised water for reconstitution e.g. MRX Laboratory Water | K5036 |
| Phosphate buffered saline (PBS) for dilution, e.g. MRX PBS Diluent | K5047 |

11 Quality control

To maintain consistent assay results, it is recommended that control plasmas are assayed at regular intervals. MRX Specialty Controls (K5021/K5022) are recommended for the assays MRX Red Free Protein S, MRX Blue Free Protein S, MRX Protein C and MRX Green vWF Antigen. Each laboratory should establish a control range to determine the allowable variation in the day-to-day performance of the test, as well as appropriate intervals for analysing controls in accordance with good laboratory practice. Recalibration is suggested, as a minimum, whenever control plasmas are not within the acceptable range and each time a new batch of reagent is used.

12 Results

For reporting of results, refer to the Instructions for Use associated with each individual assay (K5009, K5010, K5020, K5052).

13 Expected values

For expected values, refer to the Instructions for Use associated with each individual assay (K5009, K5010, K5020, K5052).

14 Limitations and interfering substances

For limitations and interfering substances, refer to the Instructions for Use associated with each individual assay (K5009, K5010, K5020, K5052).

15 Analytical performance characteristics

For analytical performance characteristics, refer to the Instructions for Use associated with each individual assay (K5009, K5010, K5020, K5052).

16 Reporting of incidents

Any serious incidents that occur in relation to this device shall be reported to Nordic Biomarker as well as the national competent authority in which the user is established.

17 Additional information

A paper copy of these Instructions for Use is available on request. Contact your local distributor.

The instrument-specific application sheet is available from your local distributor.

18 Definition of symbols

Manufacturer

Consult electronic instructions for use

nordicbiomarker.com/IFU

CE mark

Use-by date

In vitro diagnostic medical device

Temperature limit

Catalogue number

Biological risks

Batch code

Contains human blood or plasma derivatives

19 Revision history

| Version | Changes to previous version |
|---------|---------------------------------|
| 7.0 | Added translation into Italian. |

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