Instructions for Use

MRX Specialty Abnormal Control, K5022

Contents

<u>English [EN]</u>	<u>3</u>
<u>Dansk [DA]</u>	<u>6</u>
<u>Deutsch [DE]</u>	<u>9</u>
<u>Español [ES]</u>	<u>12</u>
<u>Suomi [FI]</u>	<u>15</u>
<u>Magyar [HU]</u>	<u>18</u>
<u>Italiano [IT]</u>	<u>21</u>
<u>Norsk [NO]</u>	<u>24</u>
<u>Svenska [SV]</u>	<u>27</u>

Instructions for Use [EN]

MRX Specialty Abnormal Control

REF K5022

For *In vitro* Diagnostic Use.

1 Intended use

For quality control of the assays MRX Red Free Protein S, MRX Blue Free Protein S, MRX Protein C and MRX Green vWF Antigen. Intended to be used by professional laboratory personnel using coagulation analysers.

2 Background and principle of method

For background information and principle of the method, refer to the Instructions for Use associated with each individual assay (K5009, K5010, K5020, K5052).

3 Components

MRX Specialty Abnormal Control consists of: 10 × 1 mL lyophilised citrated human plasma with additives of protein components from bovine plasma and preservative.

4 Metrological traceability

Each new lot of MRX Specialty Abnormal Control is assigned with each corresponding MRX assay against calibrators with traceability to WHO International Standards, NIBSC code: 03/228 (Free Protein S), 02/342 (Protein C) and 07/316 (vWF Antigen). Refer to the Certificate of Analysis for the lot-specific concentration of each individual analyte.

Specifications for each analyte are listed below:

Analyte	Specification
Free Protein S*	35–55%
Protein C	35–55%
vWF Antigen	20–50%

*MRX Blue Free Protein S is used for Free Protein S assignment.

5 Warnings and precautions

Wear suitable clothing for protection. Avoid contact with skin and eyes. Do not empty into drains. Waste must be disposed of in accordance with local regulations.

The control contains material of human origin. Each donor has been tested by approved methods and found negative for the presence of HBsAg and anti-HIV I & II and anti-HCV. However, as no method can offer complete assurance that infectious agents are absent, this material should be handled as potentially infectious.

The control contains Bovine Serum Albumin (< 4.0%) for stability, and bovine plasma. The animals were approved by veterinarians by ante- and post-mortem inspections. However, as no method can offer complete assurance, this material should be handled as potentially infectious.

The control contains sodium azide (less than 0.001%) to prevent microbial growth; use proper disposal procedures.

6 Preparation

- Before opening, carefully tap the vial against a surface to collect the lyophilised material at the bottom.
- Add 1.00 mL deionised water (e.g. MRX Laboratory Water, K5036). The water temperature should be 15 - 25 °C.
- Reseal the vial and let it stand for approximately 15 minutes at 15 - 25 °C.
- Gently mix by swirling or rotating until the content is completely reconstituted.

7 Storage and stability

Store at 2 - 8 °C. After reconstitution, stable for 24 hours at 2 - 25 °C in the closed original vial, provided no contamination occurs.

8 Specimen collection and preparation

For specimen collection and preparation, refer to the Instructions for Use associated with each individual assay (K5009, K5010, K5020, K5052).

9 Procedure

For description of procedure, refer to the Instructions for Use associated with each individual assay (K5009, K5010, K5020, K5052).

10 Material required but not provided

Coagulation analyser, pipettes and the following:

Reagent	REF
MRX Red Free Protein S	K5009
MRX Blue Free Protein S	K5010
MRX Protein C	K5020
MRX Green vWF Antigen	K5052
Control material	REF
MRX Specialty Normal Control	K5021
Solutions	REF
Deionised water for reconstitution e.g. MRX Laboratory Water	K5036
Phosphate buffered saline (PBS) for dilution, e.g. MRX PBS Diluent	K5047

11 Quality control

To maintain consistent assay results, it is recommended that control plasmas are assayed at regular intervals. MRX Specialty Controls (K5021/K5022) are recommended for the assays MRX Red Free Protein S, MRX Blue Free Protein S, MRX Protein C and MRX Green vWF Antigen. Each laboratory should establish a control range to determine the allowable variation in the day-to-day performance of the test, as well as appropriate intervals for analysing controls in accordance with good laboratory practice. Recalibration is suggested, as a minimum, whenever control plasmas are not within the acceptable range and each time a new batch of reagent is used.

12 Results

For reporting of results, refer to the Instructions for Use associated with each individual assay (K5009, K5010, K5020, K5052).

13 Expected values

For expected values, refer to the Instructions for Use associated with each individual assay (K5009, K5010, K5020, K5052).

14 Limitations and interfering substances

For limitations and interfering substances, refer to the Instructions for Use associated with each individual assay (K5009, K5010, K5020, K5052).

15 Analytical performance characteristics

For analytical performance characteristics, refer to the Instructions for Use associated with each individual assay (K5009, K5010, K5020, K5052).

16 Reporting of incidents

Any serious incidents that occur in relation to this device shall be reported to Nordic Biomarker as well as the national competent authority in which the user is established.

17 Additional information

A paper copy of these Instructions for Use is available on request. Contact your local distributor.

The instrument-specific application sheet is available from your local distributor.

18 Definition of symbols

Manufacturer

Use-by date

CE mark

Temperature limit

In vitro diagnostic medical device

Biological risks

Catalogue number

Contains biological material of animal origin

Batch code

Contains human blood or plasma derivatives

Consult electronic instructions for use

nordicbiomarker.com/IFU

19 Revision history

Version	Changes to previous version
6.0	Added translation into Italian. Section 18: Added symbol for biological material of animal origin. (HU)

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