Instructions for Use

MRX Protein C Calibrator, K5023

Contents

English [EN] _____ 3

Español [ES] _____ 6

Italiano [IT] _____ 9

Instructions for Use [EN]

MRX Protein C Calibrator

REF K5023

For *In vitro* Diagnostic Use.

1 Intended use

For calibration of the MRX Protein C reagent in percent. Intended to be used by professional laboratory personnel using coagulation analysers.

2 Background and principle of method

For background information and principle of the method, refer to the Instructions for Use associated with MRX Protein C (K5020).

3 Components

MRX Protein C Calibrator consists of:
4 × 1 mL lyophilised citrated human plasma with specified level of Protein C.

4 Metrological traceability

Each new lot of MRX Protein C Calibrator is assigned with MRX Protein C against in-house reference material with traceability to WHO International Standard (NIBSC code: 02/342). Refer to the Certificate of Analysis or the label for the lot-specific Protein C concentration.

Specification:

Analyte	Specification
Protein C	80 - 120%

5 Warnings and precautions

Wear suitable clothing for protection. Avoid contact with skin and eyes. Do not empty into drains. Waste must be disposed of in accordance with local regulations.

The calibrator contains material of human origin. Each donor has been tested by approved methods and found negative for the presence of HBsAg and anti-HIV I & II and anti-HCV. However, as no method can offer complete assurance that infectious agents are absent, this material should be handled as potentially infectious.

6 Preparation

- Before opening, carefully tap the vial against a surface to collect the lyophilised material at the bottom.
- Add 1.00 mL deionised water (e.g. MRX Laboratory Water, K5036). The water temperature should be 15 - 25 °C.
- Reseal the vial and let it stand for approximately 15 minutes at 15 - 25 °C.
- Gently mix by swirling or rotating until the content is completely reconstituted.

7 Storage and stability

Store at 2 - 8 °C. After reconstitution, stable for 24 hours at 2 - 25 °C in the closed original vial, provided no contamination occurs.

8 Specimen collection and preparation

For specimen collection and preparation, refer to the Instructions for Use associated MRX Protein C (K5020).

9 Procedure

For each instrument, refer to its operator's manual and to the instrument-specific application sheet.

10 Material required but not provided

Coagulation analyser capable of detection 405 nm, pipettes and the following:

Reagent	REF
MRX Protein C	K5020

Control material	REF
MRX Specialty Normal Control	K5021
MRX Specialty Abnormal Control	K5022

Solutions	REF
Deionised water for reconstitution e.g. MRX Laboratory Water	K5036
Phosphate buffered saline (PBS) for dilution, e.g. MRX PBS Diluent	K5047

11 Quality control

To maintain consistent assay results, it is recommended that control plasmas are assayed at regular intervals. MRX Specialty Controls (K5021/K5022) are recommended for MRX Protein C. Each laboratory should establish a control range to determine the allowable variation in the day-to-day performance of the test, as well as appropriate intervals for analysing controls in accordance with good laboratory practice. Recalibration is suggested, as a minimum, whenever control plasmas are not within the acceptable range and each time a new batch of reagent is used.

12 Results

The results are reported in % Protein C, where 100% corresponds to 1 IU/mL Protein C of the WHO 2nd International Standard (NIBSC code 02/342).

13 Expected values

For expected values, refer to the Instructions for Use associated with MRX Protein C (K5020).

14 Limitations and interfering substances

For limitations and interfering substances, refer to the Instructions for Use associated with MRX Protein C (K5020).

15 Analytical performance characteristics

For analytical performance characteristics, refer to the Instructions for Use associated with MRX Protein C (K5020).

16 Reporting of incidents

Any serious incidents that occur in relation to this device shall be reported to Nordic Biomarker as well as the national competent authority in which the user is established.

17 Additional Information

A paper copy of these Instructions for Use is available on request. Contact your local distributor.

The instrument-specific application sheet is available from your local distributor.

18 Definition of symbols

Manufacturer

Consult electronic instructions for use

nordicbiomarker.com/IFU

CE mark

Use-by date

In vitro diagnostic medical device

Temperature limit

Catalogue number

Biological risks

Batch code

Contains human blood or plasma derivatives

19 Revision history

Version	Changes to previous version
5.0	Added translation into Italian.

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