Instructions for Use

MRX PT Quick, K5024

Contents

| | |
|----------------------|----|
| <u>English [EN]</u> | 3 |
| <u>Español [ES]</u> | 7 |
| <u>Italiano [IT]</u> | 11 |

Instructions for Use [EN]

MRX PT Quick

REF K5024

For *In vitro* Diagnostic Use.

1 Intended use

Clotting test for quantitative determination of the prothrombin time (PT) in citrated human plasma according to the Quick method. Can be used to monitor warfarin therapy and to evaluate the function of the extrinsic and common pathway of the coagulation cascade. Intended to be used by professional laboratory personnel using coagulation analysers.

2 Background and principle of method

The MRX PT Quick reagent is used for determination of PT in patients administered vitamin K-antagonists and for detection of deficiencies within the extrinsic and common pathways of the coagulation cascade. MRX PT Quick is sensitive to abnormalities of coagulation factors II (prothrombin), V, VII, X, and fibrinogen (factor I) in citrated plasma.¹ The PT Quick method is based on stimulation of the extrinsic coagulation pathway via activation of FVII by thromboplastin in the presence of calcium. The activated complex (FVIIa-tissue factor) activates FX. Subsequently, FXa, in complex with FV, activates prothrombin to thrombin leading to the conversion of fibrinogen to fibrin, detected as a clot. In patients treated with vitamin K-antagonists, the activity of vitamin K-dependent coagulation factors will be lower and result in a prolonged PT.

3 Components

MRX PT Quick consists of lyophilised rabbit brain thromboplastin, buffer with calcium chloride, stabiliser, and preservative.

MRX PT Quick is available in the following packaging sizes:

| REF | No of vials | Volume |
|-------|-------------|--------|
| K5024 | 10 | 5 mL |

4 Warnings and precautions

Wear suitable clothing for protection. Avoid contact with skin and eyes. Do not empty into drains. Waste must be disposed of in accordance with local regulations.

MRX PT Quick contains sodium azide (less than 0.1%) to prevent microbial growth; use proper disposal procedures.

5 Preparation

- Allow the vial of MRX PT Quick to equilibrate at 15 - 25 °C for at least 10 minutes before reconstitution.
- Dissolve the content of each vial in deionised water (e.g. MRX Laboratory Water, K5036) according to:
 - MRX PT Quick K5024 in 5 mL
- Swirl gently. Keep the reconstituted reagent at 15 - 25 °C for 60 minutes and mix by inverting the vial before use. Make sure that the product is completely reconstituted before use. The reagent will dissolve into a slightly opaque colourless liquid.
- Continuous stirring or repeated inversion of the reagent is necessary during analysis.

6 Storage and stability

Store at 2 - 8 °C. After reconstitution, stable for 5 days at 2 - 25 °C in the closed original vial, provided no contamination occurs.

On-board stability: 4 days at 15 °C in open vial on ACL Top 700 CTS.

7 Specimen collection and preparation

Venous blood is collected in 3.2% sodium citrate at a ratio of 9 parts blood to 1 part anticoagulant (1:10 ratio). The ratio is critical. Trauma or stasis during blood sampling should be avoided. Blood should not be collected through a heparin lock or other heparinised line. Inverse immediately after sampling. The presence of any clots in a specimen is a cause for rejection. Centrifuge to produce platelet-poor

plasma and use for analysis. Refer to CLSI guideline H21-A5 for further instructions on specimen collection, handling and storage.²

8 Procedure

8.1 Calibration

Each lot of MRX PT Quick must be calibrated prior to determination of PT/INR in patients' plasma samples. During calibration, reagent lot and instrument specific values for normal clotting time (NCT) and International Sensitivity Index (ISI) required for INR determination (see Section 11 Results) are established.

Local NCT and ISI calibration can be performed using lyophilised calibrant plasmas,³ such as MRX PT Quick (INR) Calibrator Set.

It is also possible to determine lot and instrument specific mean normal prothrombin time (MNPT) and ISI using guidelines provided by WHO.⁴

8.2 Automated instruments

For each instrument, refer to its operator's manual and to the instrument-specific application sheet.

8.3 Semi-automated and manual instruments

- Make sure that the final volume in the cuvette is within instrument specification.
- Preincubate the reconstituted reagent at 37 °C.
- Incubate 100 µL of plasma at 37 °C for 1 minute.
- Thoroughly mix the pre-incubated MRX PT Quick reagent and add 200 µL of the reagent to the plasma sample and start recording the time.
- Record the clotting time in seconds.

9 Material required but not provided

Coagulation analyser, cuvettes, pipettes, and the following:

| Calibrator | REF |
|-------------------------------|-------|
| MRX PT Quick (INR) Calibrator | K5043 |

| Control material | REF |
|------------------------------|-------|
| MRX Routine Normal Control | K5039 |
| MRX Routine Abnormal Control | K5040 |

| Solutions | REF |
|---|-------|
| Deionised water for reconstitution e.g. MRX Laboratory Water | K5036 |

10 Quality control

To maintain consistent assay results, it is recommended that control plasmas are assayed at regular intervals. MRX Routine Controls (K5039/K5040) are recommended for MRX PT Quick. Each laboratory should establish a control range to determine the allowable variation in the day-to-day performance of the test, as well as appropriate intervals for analysing controls in accordance with good laboratory practice.

11 Results

PT Quick results are expressed in either seconds or international normalised ratio (INR). The international recommendation of reporting the PT time is INR. INR is calculated by dividing the obtained clotting time by the normal clotting time or mean normal PT (see Section 8.1 Calibration) raised to the lot specific ISI value for the reagent and measurement system according to the formula below.^{5,6}

$$\text{INR} = (\text{patient's PT} / \text{control PT})^{\text{ISI}}$$

Patient's PT= The patient's PT in seconds.

Control PT= The normal PT time (NCT or MNPT) in seconds.

Lot specific value for the reagent and instrument combination determined during calibration.

ISI = International Sensitivity Index, specific for each batch of reagent and instrument and reagent system.

Note: The ISI and NCT/MNPT values are specific for each batch of reagent and each reagent/instrument combination.

12 Expected values

INR

Expected values for non-affected population are INR 0.9 - 1.2⁷ while expected values for affected population (e.g. patients treated with vitamin K-antagonists) are INR 2.0 - 3.5.⁸

Seconds

Results from PT measurements in plasma samples from 95 healthy blood donors using MRX PT Quick are presented below. The analysis was performed using an ACL Top 700 CTS instrument.

| N | Mean PT | Mean ± 2 SD |
|------------|---------|---------------|
| 95 samples | 12.6 s | 11.1 - 14.1 s |

Due to inter-laboratory variation, each laboratory should set its own reference interval.

13 Limitations and interfering substances

Results may be affected by insufficient blood sampling with shifted ratio of sodium citrate to patient plasma or by interfering substances such as heparin and EDTA.

MRX PT Quick is affected by direct/novel oral anticoagulants (DOACs/NOACs). Therefore, analysis of DOAC containing plasma samples is not recommended using MRX PT Quick.

MRX PT Quick is insensitive to the following substances on ACL Top 700 CTS:

| Interfering substance | Tolerance |
|------------------------|-----------------|
| Bilirubin | Up to 40 mg/dL |
| Haemoglobin | Up to 500 mg/dL |
| Triglycerides | Up to 500 mg/dL |
| Unfractionated heparin | Up to 200 U/dL |

14 Analytical performance characteristics

The following performance data was obtained with an ACL Top 700 CTS instrument. Performance will depend on the instrument used.

When compared to other PT Quick reagents, MRX PT Quick correlates as follows:

$$y \text{ (MRX PT Quick on ACL Top 700 CTS)} = 1.4 \times (\text{STA Neoplastine CI Plus on ACL Top 700 CTS}) - 0.6; r^2 = 0.97.$$

$$y \text{ (MRX PT Quick on ACL Top 700 CTS)} = 1.2 \times (\text{RecombiPlasTin 2 G on ACL Top 700 CTS}) - 0.2; r^2 = 0.97.$$

Precision:

| Sample | Mean PT | Repeatability CV |
|---------|---------|------------------|
| Level 1 | 13.7 s | 0.7% |
| Level 2 | 29.4 s | 2.0% |

15 Reporting of incidents

Any serious incidents that occur in relation to this device shall be reported to Nordic Biomarker as well as the national competent authority in which the user is established.

16 Additional information

A paper copy of these Instructions for Use is available on request. Contact your local distributor.

The instrument-specific application sheet is available from your local distributor.

17 References

- QUICK, Armand J. The thromboplastin reagent for the determination of prothrombin. *Science*, 1940, 92.2379: 113-114.
- CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline – Fifth Edition. CLSI document H21-A5. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
- VAN DEN BESSELAAR, A. M. H. P., et al. Guidelines on preparation, certification, and use of certified plasmas for ISI calibration and INR determination. *Journal of Thrombosis and Haemostasis*, 2004, 2.11: 1946-1953.
- WHO EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION. MEETING; WORLD HEALTH ORGANIZATION; WORLD HEALTH ORGANIZATION. EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION. *Who Expert Committee on Biological Standardization: Sixty-second Report*. World Health Organization, 2013.
- LOELINGER, E. A. ICSH/ICTH recommendation for reporting prothrombin time in oral anticoagulant control. *Thromb Haemostas*, 1985, 53: 155-156.
- WHO EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION, et al. *WHO Expert Committee on Biological Standardization: forty-eighth report*. World Health Organization, 1999.
- ZIERK, Jakob, et al. Data mining of reference intervals for coagulation screening tests in adult patients. *Clinica Chimica Acta*, 2019, 499: 108-114.
- HIRSH, Jack, et al. Oral anticoagulants: mechanism of action, clinical effectiveness, and optimal therapeutic range. *Chest*, 1998, 114.5: 445S-469S.

18 Definition of symbols

Manufacturer

Consult electronic instructions for use

nordicbiomarker.com/IFU

CE mark

Use-by date

In vitro diagnostic medical device

Temperature limit

Catalogue number

Biological risks

Batch code

Contains biological material of animal origin

19 Revision history

| Version | Changes to previous version |
|---------|---------------------------------|
| 6.0 | Added translation into Italian. |

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17 Referencias

- QUICK, Armand J. The thromboplastin reagent for the determination of prothrombin. *Science*, 1940, 92.2379: 113-114.
- CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline – Fifth Edition. CLSI document H21-A5. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
- VAN DEN BESSELAAR, A. M. H. P., et al. Guidelines on preparation, certification, and use of certified plasmas for ISI calibration and INR determination. *Journal of Thrombosis and Haemostasis*, 2004, 2.11: 1946-1953.
- WHO EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION. MEETING; WORLD HEALTH ORGANIZATION; WORLD HEALTH ORGANIZATION. EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION. *Who Expert Committee on Biological Standardization: Sixty-second Report*. World Health Organization, 2013.
- LOELINGER, E. A. ICSH/ICTH recommendation for reporting prothrombin time in oral anticoagulant control. *Thromb Haemostas*, 1985, 53: 155-156.

6. WHO EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION, et al. *WHO Expert Committee on Biological Standardization: forty-eighth report*. World Health Organization, 1999.
7. ZIERK, Jakob, et al. Data mining of reference intervals for coagulation screening tests in adult patients. *Clinica Chimica Acta*, 2019, 499: 108-114.
8. HIRSH, Jack, et al. Oral anticoagulants: mechanism of action, clinical effectiveness, and optimal therapeutic range. *Chest*, 1998, 114.5: 445S-469S.

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17 Bibliografia

1. QUICK, Armand J. The thromboplastin reagent for the determination of prothrombin. *Science*, 1940, 92.2379: 113-114.
2. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline – Fifth Edition. CLSI document H21-A5. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
3. VAN DEN BESSELAAR, A. M. H. P., et al. Guidelines on preparation, certification, and use of certified plasmas for ISI calibration and INR determination. *Journal of Thrombosis and Haemostasis*, 2004, 2.11: 1946-1953.
4. WHO EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION. MEETING; WORLD HEALTH ORGANIZATION; WORLD HEALTH ORGANIZATION. EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION. *Who Expert Committee on Biological Standardization: Sixty-second Report*. World Health Organization, 2013.
5. LOELINGER, E. A. ICSH/ICTH recommendation for reporting prothrombin time in oral anticoagulant control. *Thromb Haemostas*, 1985, 53: 155-156.
6. WHO EXPERT COMMITTEE ON BIOLOGICAL STANDARDIZATION, et al. *WHO Expert Committee on Biological Standardization: forty-eighth report*. World Health Organization, 1999.

7. ZIERK, Jakob, et al. Data mining of reference intervals for coagulation screening tests in adult patients. *Clinica Chimica Acta*, 2019, 499: 108-114.
8. HIRSH, Jack, et al. Oral anticoagulants: mechanism of action, clinical effectiveness, and optimal therapeutic range. *Chest*, 1998, 114.5: 445S-469S.

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