# Instructions for Use

MRX APTT, K5029, K5030

## Contents

English [EN] \_\_\_\_\_ 3

Español [ES] \_\_\_\_\_ 7

Italiano [IT] \_\_\_\_\_ 11

Lietuvių [LT] \_\_\_\_\_ 15

# Instructions for Use [EN]

## MRX APTT

**REF** K5029, K5030

For *In vitro* Diagnostic Use.

### 1 Intended use

Clotting test for quantitative determination of the Activated Partial Thromboplastin Time (APTT) in citrated human plasma. Can be used to evaluate the function of the intrinsic and common pathway of the coagulation cascade. Intended to be used by professional laboratory personnel using coagulation analysers.

### 2 Background and principle of method

In the APTT test, a contact activator is used to stimulate the generation of Factor XIa by providing a surface for the interaction of FXI, high molecular weight kininogen, kallikrein and Factor XIIa.<sup>1</sup> The contact activation is allowed to proceed at 37 °C for a specific period. Calcium chloride is then added to trigger further reactions and the time required for formation of a fibrin clot is determined. Phospholipids are required to form complexes which activate Factor X and prothrombin.

The APTT reagent contains phospholipids and silica to ensure a highly consistent and stable product.<sup>2</sup> Furthermore, the APTT reagent is lupus anticoagulant insensitive. Lupus anticoagulant insensitive reagents yield more reliable factor assay results than reagents which are sensitive to lupus inhibitors.<sup>3,4</sup> Lupus inhibitors occur occasionally among test plasmas.

Prolonged clotting times may be observed in the following situations: deficiency of Factor XII, XI, X, IX; VIII; V, II or fibrinogen, in liver diseases, vitamin K deficiency, presence of heparin, lupus anticoagulant or other inhibitors, which affect the intrinsic pathway.<sup>5</sup> A shortened APTT appears to correlate with an increased risk of venous thrombosis.<sup>6</sup>

### 3 Components

MRX APTT consists of colloidal silicate with phospholipids, buffer, lactose from bovine milk, and preservatives.

MRX APTT is available in the following packaging sizes:

REF	No of vials	Volume
K5029	10	5 mL
K5030	10	10 mL

### 4 Warnings and precautions

Wear suitable clothing for protection. Avoid contact with skin and eyes. Do not empty into drains. Waste must be disposed of in accordance with local regulations.

The reagent contains sodium azide (less than 0.1%) to prevent microbial growth; use proper disposal procedures.

### 5 Preparation

Mix by inverting the vial before use.

### 6 Storage and stability

Store at 2 - 8 °C. After opening, stable for 30 days at 2 - 8 °C in the closed original vial, provided no contamination occurs.

### 7 Specimen collection and preparation

Venous blood is collected in 3.2% sodium citrate at a ratio of 9 parts blood to 1 part anticoagulant (1:10 ratio). The ratio is critical. Trauma or stasis during blood sampling should be avoided. Blood should not be collected through a heparin lock or other heparinised line. Inverse immediately after sampling. The presence of any clots in a specimen is a cause for rejection. Centrifuge to produce platelet-poor plasma and use for analysis. Refer to CLSI guideline H21-A5 for further instructions on specimen collection, handling and storage.<sup>7</sup>

## 8 Procedure

### 8.1 Automated instruments

For each instrument, refer to its operator's manual and to the instrument-specific application sheet.

### 8.2 Semi-automated and manual instruments

- Make sure that the final volume in the cuvette is within instrument specification.
- Preincubate 25 mM calcium chloride (e.g., MRX Calcium Chloride, K5049/K5050/K5051) at 37 °C for at least 10 minutes.
- Pipette 100 µL of plasma sample into a test cuvette. Incubate at 37 °C for 1 - 2 minutes.
- Add 100 µL MRX APTT reagent to the cuvette containing the plasma sample. Incubate the mixture at 37 °C for 5 minutes.
- Add 100 µL of preheated calcium chloride and start recording the time.
- Record the clotting time in seconds.

## 9 Material required but not provided

Coagulation analyser, cuvettes, pipettes, and the following:

Control material	REF
MRX Routine Normal Control	K5039
MRX Routine Abnormal Control	K5040

Solutions	REF
Calcium Chloride (25 mM), e.g. MRX	K5049
Calcium Chloride	K5050
	K5051
Deionised water for reconstitution e.g. MRX Laboratory Water	K5036

## 10 Quality control

To maintain consistent assay results, it is recommended that control plasmas are assayed at regular intervals. MRX Routine Controls (K5039/K5040) are recommended for MRX APTT. Each laboratory should establish a control range to determine the allowable variation in the day-to-day performance of the test, as well as appropriate intervals for analysing controls in accordance with good laboratory practice.

## 11 Results

The results are reported in seconds.

## 12 Expected values

Due to inter-laboratory variation, each laboratory should set its own reference interval.

## 13 Limitations and interfering substances

Delay in testing and difficulty in specimen collection may result in falsely prolonged APTT results. Hence, a strict control of plasma preparation and handling is of importance. Inadvertent pre-activation of coagulation factors may result in shortening of the APTT. The impact of coagulation factor deficiency on the APTT may be compensated for by pre-activation and by elevated levels of other coagulation factors.

APTT results may be affected by drugs targeting coagulation factors in the intrinsic and common coagulation pathways.

MRX APTT is insensitive to the following substances on ACL9000:

Interfering substance	Tolerance
Bilirubin	40 mg/dL
Haemoglobin	1000 mg/dL

Heparin sensitivity:

Heparin U/mL	APTT clotting time in seconds
0	29.9
0.2	70.0
0.4	174

Heparin sensitivity was evaluated using PNP spiked with unfractionated heparin. All tests were run with 5 minutes preincubation on ACL300R. Due to many variables (e.g. different sources of heparin) which may affect the clotting times, each laboratory should establish its own instrument-specific level where heparin is detected in citrated plasma.

Factor sensitivity:

Factor (%)	APTT-Clotting Time (in seconds)		
	Expected values 26.7 - 34.9 sec.		
	FVIII	FIX	FXI
<1	85.7	70.9	92.4
10	45.9	44.4	50.7
40	34.2	33.9	34.4
100	28.8	28.8	28.8

Factor sensitivity testing was carried out by mixing varying volumes of pooled normal plasma (PNP, 100% factor levels) with individual intrinsic factor deficient plasma to achieve 10% and 40% factor levels. Results with PNP and factor deficient plasmas alone are also shown in the table above. All tests were run with 5 minutes preincubation on ACL300R. Expected values shown are mean  $\pm$  2 SD. These values should be used as guidelines only. Each laboratory should establish factor sensitivity using their own instruments and techniques.

## 14 Analytical performance characteristics

The following performance data was obtained with an ACL 9000 instrument. Performance will depend on the instrument used.

Precision:

Sample	Mean APTT	Repeatability CV
Level 1	35.6	0.9%
Level 2	71.3	1.2%

## 15 Reporting of incidents

Any serious incidents that occur in relation to this device shall be reported to Nordic Biomarker as well as the national competent authority in which the user is established.

## 16 Additional information

A paper copy of these Instructions for Use is available on request. Contact your local distributor.

The instrument-specific application sheet is available from your local distributor.

## 17 References

1. PROCTOR, Robert R.; RAPAPORT, Samuel I. The partial thromboplastin time with kaolin: a simple screening test for first stage plasma clotting factor deficiencies. *American Journal of Clinical Pathology*, 1961, 36.3: 212-219.
2. STEVENSON, K. J., et al. The reliability of activated partial thromboplastin time methods and the relationship to lipid composition and ultrastructure. *Thrombosis and haemostasis*, 1986, 55.02: 250-258.
3. BRANDT, J. T., et al. Effect of lupus anticoagulants on the activated partial thromboplastin time. Results of the College of American Pathologists survey program. *Archives of pathology & laboratory medicine*, 1991, 115.2: 109-114.
4. DENIS-MAGDELAINE, A.; FLAHAULT, Antoine; VERDY, Elisabeth. Sensitivity of sixteen APTT reagents for the presence of lupus anticoagulants. *Pathophysiology of Haemostasis and Thrombosis*, 1995, 25.3: 98-105.
5. HATHAWAY, William E., et al. Activated partial thromboplastin time and minor coagulopathies. *American journal of clinical pathology*, 1979, 71.1: 22-25.
6. TRIPODI, Armando, et al. A shortened activated partial thromboplastin time is associated with the risk of venous thromboembolism. *Blood*, 2004, 104.12: 3631-3634.
7. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline – Fifth Edition. CLSI document H21-A5. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.

## 18 Definition of symbols

Manufacturer

Consult electronic instructions for use

[nordicbiomarker.com/IFU](http://nordicbiomarker.com/IFU)

CE mark

Use-by date

In vitro diagnostic medical device

Temperature limit

Catalogue number

Contains biological material of animal origin

Batch code

## 19 Revision history

Version	Changes to previous version
4.0	Added translation into Spanish.

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