# Instructions for Use

MRX Antithrombin, K5033

# Contents

|               |    |
|---------------|----|
| English [EN]  | 3  |
| Dansk [DA]    | 6  |
| Deutsch [DE]  | 10 |
| Español [ES]  | 14 |
| Suomi [FI]    | 18 |
| Italiano [IT] | 22 |
| Norsk [NO]    | 26 |
| Svenska [SV]  | 30 |

# Instructions for Use [EN]

## MRX Antithrombin

REF

 K5033

For *In vitro* Diagnostic Use.

### 1 Intended use

Chromogenic test for quantitative determination of antithrombin (AT) activity in citrated human plasma as an aid in diagnosis of antithrombin deficiency. Intended to be used by professional laboratory personnel using coagulation analysers measuring at 405 nm.

### 2 Background and principle of method

Antithrombin (AT) is the primary physiological inhibitor of thrombin and factor Xa (FXa) in plasma and thereby effectively regulates blood coagulation. Unfractionated (UF) and low molecular weight (LMW) heparin greatly enhances AT activity. Hereditary or acquired AT deficiency is an important risk factor for venous thromboembolic disorders.<sup>1,2</sup> MRX Antithrombin is not affected by Heparin Cofactor II activity.

Antithrombin activity is determined in a two-stage chromogenic assay as follows:

- An excess of FXa is added to citrated plasma in the presence of heparin.
- The residual FXa activity is determined from the hydrolysis of a chromogenic FXa-substrate, which results in the release of free p-nitroaniline, pNA. The AT activity is inversely related to the amount of released pNA and is expressed as % AT activity derived from a calibration curve.

### 3 Components

MRX Antithrombin consists of:

- Bovine Factor Xa: 3 × 6 mL of bovine FXa, 10 nkat/mL, in buffer containing heparin, stabilisers and preservatives.
- Factor Xa Substrate: 3 × 4.5 mL of chromogenic FXa substrate in water medium containing surfactant and preservatives.

### 4 Warnings and precautions

Wear suitable clothing for protection. Avoid contact with skin and eyes. Do not empty into drains. Waste must be disposed of in accordance with local regulations.

Bovine Factor Xa contains Bovine Serum Albumin, for stability and bovine FXa. The animals were approved by veterinarians by ante- and post-mortem inspections. However, as no method can offer complete assurance, this material should be handled as potentially infectious.

### 5 Preparation

- Bovine Factor Xa: Ready to use. Mix by inverting the vial before use.
- Factor Xa Substrate: Ready to use. Mix by inverting the vial before use.

### 6 Storage and stability

- Bovine Factor Xa: Store at 2 - 8 °C. After opening, stable for 1 month at 2 - 8 °C and for 48 hours at 12 - 15 °C in the closed original vial, provided no contamination occurs.
- Factor Xa Substrate: Store at 2 - 8 °C. After opening, stable for 1 month at 2 - 8 °C and for 48 hours at 12 - 15 °C in the closed original vial, provided no contamination occurs.  
Note: Factor Xa Substrate is light-sensitive and should be stored in the dark.

### 7 Specimen collection and preparation

Venous blood is collected in 3.2% sodium citrate at a ratio of 9 parts blood to 1 part anticoagulant (1:10 ratio). The ratio is critical. Trauma or stasis during blood sampling should be avoided. Blood should not be collected through a heparin lock or other heparinised line. Inverse immediately after sampling. The presence of any clots in a specimen is a cause for rejection. Centrifuge to produce platelet-poor plasma and use for analysis. Refer to CLSI guideline H21-A5 for further instructions on specimen collection, handling and storage.<sup>3</sup>

## 8 Procedure

For each instrument, refer to its operator's manual and to the instrument-specific application sheet.

## 9 Material required but not provided

Coagulation analyser capable of detection 405 nm, cuvettes, pipettes, and the following:

|                                                                |       |
|----------------------------------------------------------------|-------|
| Calibrator                                                     | REF   |
| MRX Antithrombin Calibrator                                    | K5041 |
| Control material                                               | REF   |
| MRX Routine Normal Control                                     | K5039 |
| MRX Routine Abnormal Control                                   | K5040 |
| Solutions                                                      | REF   |
| Saline solution for dilution, e.g. MRX Sodium Chloride Diluent | K5046 |
| Deionised water for reconstitution e.g. MRX Laboratory Water   | K5036 |

## 10 Quality control

To maintain consistent assay results, it is recommended that control plasmas are assayed at regular intervals. MRX Routine Controls (K5039/K5040) are recommended for MRX Antithrombin. Each laboratory should establish a control range to determine the allowable variation in the day-to-day performance of the test, as well as appropriate intervals for analysing controls in accordance with good laboratory practice. Recalibration is suggested, as a minimum, whenever control plasmas are not within the acceptable range and each time a new batch of reagent is used.

## 11 Results

The results are reported in IU/mL. 100% AT activity is defined as 1 IU AT function in 1 mL plasma of the WHO 3<sup>rd</sup> International Standard (NIBSC code 08/258).

Samples that are reported above the measuring range should be manually diluted and re-analysed. No result outside the measuring range should be used in forming a diagnosis or for patient management.

## 12 Expected values

Expected values for antithrombin activity are typically in the 0.84 - 1.17 IU/mL range.<sup>4</sup>

Results from antithrombin activity measurements in plasma samples from 50 healthy blood donors using MRX Antithrombin are presented below. The analysis was performed using an ACL Top 700 CTS instrument.

| N          | Mean AT | Mean $\pm$ 2 SD   |
|------------|---------|-------------------|
| 50 samples | 1.00    | 0.80 - 1.20 IU/mL |

Due to inter-laboratory variation, each laboratory should set its own reference interval.

## 13 Limitations and interfering substances

MRX Antithrombin is insensitive to the following substances on ACL 9000:

| Interfering substance        | Tolerance       |
|------------------------------|-----------------|
| Bilirubin                    | Up to 40 mg/dL  |
| Haemoglobin                  | Up to 150 mg/dL |
| Triglycerides                | Up to 500 mg/dL |
| Unfractionated heparin       | Up to 330 U/dL  |
| Low molecular weight heparin | Up to 330 U/dL  |

## 14 Analytical performance characteristics

The following performance data was obtained with an ACL 9000 and ACL Top 700 instrument. Performance will depend on the instrument used.

MRX Antithrombin has a measuring range of 0.2 - 1.4 IU/mL. When compared to another chromogenic assay, MRX Antithrombin correlates as follows:

$y$  (MRX Antithrombin on ACL Top 700 CTS) =  $1.02 \times$  (IL HemosIL AT on ACL Top) + 0.7;  $r^2 = 0.97$ .

Precision:

| Sample  | Mean AT    | Repeatability CV |
|---------|------------|------------------|
| Level 1 | 0.45 IU/mL | 2.6%             |
| Level 2 | 0.90 IU/mL | 2.1%             |

## 15 Reporting of incidents

Any serious incidents that occur in relation to this device shall be reported to Nordic Biomarker as well as the national competent authority in which the user is established.

## 16 Additional information

A paper copy of these Instructions for Use is available on request. Contact your local distributor.

The instrument-specific application sheet is available from your local distributor.

## 17 References

- QUINSEY, Noelene S., et al. Antithrombin: in control of coagulation. *The international journal of biochemistry & cell biology*, 2004, 36.3: 386-389.
- DAHLBÄCK, Björn. Blood coagulation. *The Lancet*, 2000, 355.9215: 1627-1632.
- CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline – Fifth Edition. CLSI document H21-A5. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
- KITCHEN, Steven., et al. Appendix 1 Normal ranges. *Practical Hemostasis and Thrombosis* 209-214 (2005)

## 18 Definition of symbols

Manufacturer

Consult electronic instructions for use

[nordicbiomarker.com/IFU](http://nordicbiomarker.com/IFU)

CE mark

Use-by date

In vitro diagnostic medical device

Temperature limit

Catalogue number

Biological risks

Batch code

Contains biological material of animal origin

## 19 Revision history

| Version | Changes to previous version                                                                |
|---------|--------------------------------------------------------------------------------------------|
| 7.0     | Section 3: Update of kit composition due to changes in number of vials and filling volume. |

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## 17 Referencer

1. QUINSEY, Noelene S., et al. Antithrombin: in control of coagulation. *The international journal of biochemistry & cell biology*, 2004, 36.3: 386-389.
2. DAHLBÄCK, Björn. Blood coagulation. *The Lancet*, 2000, 355.9215: 1627-1632.
3. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline – Fifth Edition. CLSI document H21-A5. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
4. KITCHEN, Steven., et al. Appendix 1 Normal ranges. *Practical Hemostasis and Thrombosis* 209-214 (2005)

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## 17 Referencias

1. QUINSEY, Noeline S., et al. Antithrombin: in control of coagulation. *The international journal of biochemistry & cell biology*, 2004, 36.3: 386-389.
2. DAHLBÄCK, Björn. Blood coagulation. *The Lancet*, 2000, 355.9215: 1627-1632.
3. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline – Fifth Edition. CLSI document H21-A5. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
4. KITCHEN, Steven., et al. Appendix 1 Normal ranges. *Practical Hemostasis and Thrombosis* 209-214 (2005)

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